

Hohenstein

Quality Label "Fit"

Objective

With the help of the independent product endorsement "Best suited to the material and functional fit", you can guarantee a consistent fit for your customers – right across all your collections.

Description

- The set of sizes will be picked up on-site by our independent inspectors
- The fit is tested both in new condition and after washing as described on the care label
- The assessment is made by the department fitting, development and pattern, using fitting models whose body measurements corresponds to the body-measurements table from the customer or "Size Germany" survey.
- If the set of sizes meets all requirements, the product can be endorsed with the Hohenstein quality label at the POS



Requirements for test samples

- A product that reaches the highest assessment level both before and after washing can be awarded the Hohenstein quality label. This assures that the product satisfies the following requirements:
- The garment must fit with regard to all length and width measurements
- It must allow freedom of movement as appropriate for its purpose and function
- The fit and appearance must not be impaired by defects in workmanship
- Children's clothing must conform with DIN EN 14 682

The scope of service

- Label printing template valid for 1 year
- Test report incl. fitting score
- If required, a clarifying meeting with try-on session

Customer benefit

- The neutral assurance of product quality by an independent institution can be shown to the end user at the POS in the form of a hang-tag
- Quality assessment in new condition and after washing
- Effective prevention of complaints at the POS

What we need from you

- A set of sizes in the original material incl. care instructions
- To extend the label for another year, our inspector picks up the latest set of sizes for re-testing before the expiry of the period of validity

Time required

- Maximum 8 working days; express service: 5 working days